
IOLMaster

(Software Version 3.xx)

Gebrauchsanweisung
User's Manual
Mode d'emploi
Instrucciones para el uso

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


This diagram illustrates the assembly of the ECG machine. The components are labeled as follows:

- 1**: The main base unit, which includes a circular platform with a central vertical rod and a control knob.
- 2**: The monitor screen, which is mounted on a swivel arm extending from the base.
- 3**: The ECG lead wires, which are connected to the base and the monitor.
- 4**: The ECG lead wires, which are connected to the base and the monitor.
- 5**: The power supply unit, which is connected to the base and the monitor.
- 6**: The ECG lead wires, which are connected to the base and the monitor.
- 7**: The ECG lead wires, which are connected to the base and the monitor.

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UMDNS-Nr.: 18-014

Knowledge of this manual is required for the operation of the instrument. Therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning the safe operation of the instrument.

The specifications are subject to change; the manual is not covered by an update service.

- © Unless expressly authorized, forwarding and duplication of this document, and the utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

Windows 95® and Windows 98® are trademarks of Microsoft Corporation, Inc.

IOLMaster® is a trademark of Carl Zeiss.

SRK® is a trademark of CTI (Computational Technology Inc.).

In general, other brand names as well as names of software and hardware products are subject to trademark or patent protection. Quoting of product names is for information only and does not represent any trademark misuse.

Copyright	1
Safety notes	4
General	4
Standards and regulations	4
Notes on installation and use	5
Safe operation	6
Electrical safety	6
Light emission from the instrument	6
Operational requirements	6
Proper instrument use	7
Intended use	7
Disposal	7
Instrument description	8
Controls	8
Optionally available accessories	9
Start-up	10
Functional description	11
Operation	12
Program control	12
General remarks	12
Operation via touchpad and keyboard	13
Operation via optional keyboard and mouse	13
Screen layout	14
Overview of buttons and shortcut keys	15
Summary of result displays in Axial Length Mode	16
Menu overview	17
User Database	18
Data backup	20
Setup submenu	23
Preparing measurements	25
Switching on	25
Entry of patient data	25
Aligning the instrument to the patient	26
Axial length measurement [ALM]	27
Measurement of corneal curvature [KER]	30
Multiple keratometer measurements	31
Anterior chamber depth measurement [ACD]	32
Determination of [White-To-White] (optional)	34
Measurement of the other eye	35
Result printout	36
IOL calculation	36
Filling in the IOL database	37

IOL calculation after refractive surgery of the cornea (optional)	38
IOL calculation	40
Optimization of lens constants	42
New patient	49
Using the patient database	49
Recalling previous measurement results	50
Erasing a patient	51
Renaming a patient	51
Data export (optional)	52
Switching off the IOLMaster	53
Evaluation of ALM results	53
Example signal curves of axial length measurements	53
Valid signal curves	53
Non valid signal curves	54
Zooming the graph display	55
Post-run editing of axial length measurements	56
SNR categories	56
Shifting the measuring cursor	57
Interpretation of axial length measurements	60
Signals from inner limiting membrane (ILM)	62
Signals from the choroid	63
Tips for keratometer measurement	65
How to adjust the measuring marks	65
Measurement errors	66
Maladjustment	66
Other findings	68
Tips for anterior chamber depth measurement	71
How to adjust the instrument	71
Measurement errors	73
Maladjustment	75
Pathological findings	78
Tips for WTW measurement (optional)	80
How to adjust the instrument	80
Care and maintenance	81
What to do in case of problems or system errors?	81
Checking the measurement functions	81
Remedy of printers trouble	83
Instrument care	84
Safety inspections	85
Technical Data	86
Manufacturer's Declaration	90
Medical Product Act	90

General

This instrument has been developed and tested in accordance with national and international regulations. A high degree of instrument safety has been ensured.

Observe all safety notes and information in this manual and on device labels. These notes are marked by the following pictograms.

Caution!

Danger to the instrument!

Disconnect the power cable before opening the instrument!

Instrument type B according to EN 60601-1 (IEC 60601-1).

Caution:

Proper operation of the instrument is imperative for its safe function. Therefore, please read these operating instructions thoroughly before using the instrument.

You may obtain further information from our service organization or authorized representatives.

Standards and regulations

- ☐ Quality Assurance of Carl Zeiss Meditec meets the requirements.
- ☐ The device complies with the EC Medical Device Directive 93/42/EEC and their national equivalent in form of the German Medical Product Act (MPA) (--> *Manufacturer's Declaration on p. 90*).

Notes on installation and use

- ☐ Do not operate the instrument
 - in explosion-risk areas
 - in the presence of flammable anesthetics and volatile solvents such as alcohol, benzine or similar chemicals.
- ☐ Do not store or use this instrument in damp rooms. Avoid dripping or splashing fluids near the instrument.
- ☐ This instrument may only be modified and repaired by service technicians employed with or authorized by Carl Zeiss Meditec.
- ☐ The manufacturer is not liable for damages caused by unauthorized persons tampering with the instrument; such tampering will also void any rights to claim under warranty.
- ☐ Use this instrument only with accessory equipment delivered by Carl Zeiss Meditec.
- ☐ The instrument may be operated only by instructed and trained personnel.
- ☐ Keep the users manuals always at hand.
- ☐ The instrument may only be opened by service technicians employed with or authorized by Carl Zeiss Meditec.
- ☐ Use printers recommended by Carl Zeiss Meditec only.
- ☐ Place the printer at a distance of at least 1.5 m from the patient's seat at the instrument.
- ☐ If you use a printer of Protection Class II (without protective earth terminal), you must make sure that a power isolation transformer is included in the power line of the printer.
- ☐ If you use a printer of Protection Class I (with protective earth terminal), you must make sure that a power isolation transformer is included in the power line of the printer or it is connected to a separate stationary power outlet of the room's electric installation.
- ☐ You can obtain the required isolation transformer via our sales organization.
- ☐ Do not use a cellular telephone and other devices not complying with the requirements of EMC Class B because signals transmitted may cause the equipment to malfunction. The effects of radio signals on medical instruments depend on different factors and are, therefore, not foreseeable. To avoid electromagnetic disturbances, the instrument has to be installed and set into operation in accordance with the Technical Manual and using exclusively the components delivered by Carl Zeiss Meditec.
- ☐ Additional movable multiple sockets or extension lines must not be connected to the system.
- ☐ For calculating the lens power with the IOL Master, do not use the constants of the manufacturer, but constants optimized for the optical biometry only (see also p. 42 and p. 60).

Safe operation

Electrical safety

- ☐ The built-in power supply unit is short-circuit-proof and does not contain externally accessible fuses.
- ☐ Electrical hazards to patients and operators are eliminated if the instrument is properly used.
- ☐ The instrument may be opened only by persons authorized by the manufacturer.

Light emission from the instrument

- ☐ The instrument emits light in the visible and infrared spectral regions. It keeps the limit values as specified by EN ISO 15004 "Ophthalmic Instruments Fundamental requirements and test methods".
- ☐ The limit values as specified for Class 1 laser devices by EN 60825-1 are kept if the instrument is operated as intended (limiting 20 axial length measurements per eye and per day).
- ☐ A maximum of 20 axial length measurements per eye and per day are permitted with the IOLMaster. With proper operation of the IOLMaster this value cannot be exceeded.

Class 1 Laser Device

Operational requirements

Please take care that the following operational requirements are met when using the IOLMaster:

- ☐ Use the power cable supplied with the instrument. If the instrument is mounted to a Carl Zeiss Meditec instrument table, it is powered through this table.
- ☐ The instrument is plugged into a power outlet that has a properly connected protective earth connection.
- ☐ All cables and plugs are in perfect working condition and secured where provided.
- ☐ If the protection is reduced or electrical cables are damaged, discontinue use and take care that the instrument is not operated inadvertently. Then, call Carl Zeiss Meditec service.
- ☐ Do not cover any ventilation grids or slots.

Proper instrument use

- ☐ Always enter patient data first (last name, first name, and date of birth).
- ☐ Before measuring the anterior chamber depth, inform the patient that the light will flicker during the measurement.
- ☐ Instantly turn off the power switch of the instrument and disconnect the power cable if uncertain problems arise.
- ☐ If the instrument does not respond to your commands, turn off the power switch and then restart.
- ☐ To exit the program when finished with use click on the button of the toolbar.
Turn instrument off via the power switch.

Intended use

The IOLMaster is to be used only for the measurement of axial length, corneal curvature, anterior chamber depth and optionally for "White-To-White" determinations of the human eye, as well as for the calculation of the required intraocular lens.

The responsibility for any other use than the one intended exclusively lies with the user.

The instrument may only be used in combination with accessory units approved by Carl Zeiss Meditec. Please consult Carl Zeiss Meditec or an authorized service representative regarding use of other accessories.

Disposal

The internal control computer of the device contains electric components. At the end of its lifetime, it has to be disposed of in compliance with relevant regulations.

Controls


Fig. 1 Controls of the IOLMaster

- 1 Joystick with release button**
for adjusting the instrument horizontally (X, Y) and vertically (Z, by turning)
- 2 Display**
- 3 Reds marks**
indicates patient's eye level for optical alignment
- 4 Forehead rest**
- 5 Connector board**
- 6 Display connector**
- 7 Device control connector**
- 8 Mouse connector**
- 9 Keyboard connector**
- 10 Left mouse button**
- 11 Right mouse button**
- 12 Touchpad for cursor movement**
- 13 Chin rest**
- 14 Holding pins for paper pads**
also used for test eye alignment
- 15 Patient alignment knob**
- 16 3,5" disk drive**
for technical service only
- 17 Identification label**
- 18 Instrument lock knob**


Optionally available accessories

- ☐ Footswitch
- ☐ Printer for result printout
- ☐ Keyboard support
- ☐ Instrument table IT 3L
- ☐ Software, Option A
- ☐ Software, Option B
- ☐ Power isolation transformer

Start-up

In general, Carl Zeiss Meditec or other authorized personnel will perform the following operations.

- ☐ Remove and unpack box containing accessories.
- ☐ Remove the instrument from the box with care. Only lift or carry the instrument via the instrument base.
- ☐ Removing transportation locks:
 - Loosen instrument lock knob (**18**, Fig. 1).
 - Turn joystick clockwise (one turn) to move the instrument upward and pull out the red plate underneath the base axis (patient-side).
 - Remove red pads from the gearboxes of the instrument base.
- ☐ Electrical connection (see Fig. 2):
 - Connect display and device cables and secure them.
 - Connect keyboard cable.
 - Optional: Connect and secure printer (to LPT 1), mouse and interface cable (to COM 1), and monitor cable (to VGA).
 - Connect power cable.

- | | |
|-----------------------------|------------------------------------|
| 1 Power switch | 6 Device connector |
| 2 Power plug | 7 Display connector |
| 3 Monitor connector | 8 Printer connector (LPT 1) |
| 4 Mouse connector | 9 Serial port (COM 1) * |
| 5 Keyboard connector | |

Fig. 2 Connector board of IOLMaster

Use the power switch to turn on the instrument. After a successful quick internal test you can start operating the instrument. Some operating parameters are factory set and may be changed via the Setup Menu.

* If connecting an external PC or an external monitor, the owner/operator must ensure to meet the safety requirements as per IEC 60601-1-1 (medical electrical systems)!

Functional description

The IOLMaster is a combined biometry instrument which measures parameters of the human eye needed for intraocular lens calculation. It measures quickly and precisely the following eye parameters:

- ☐ Axial length (distance between tear film and retina; see also page 60)
- ☐ Corneal curvature
- ☐ Anterior chamber depth
- ☐ "White-To-White" (optional)

All measurements are non-contact measurements providing excellent patient comfort.

The axial length measurement is based on a patented interference optical method known as Partial Coherence Interferometry (PCI). The displayed results of the axial length measurements are comparable to ultrasonic immersion measurements of axial length via the use of an internal, statistically verified calculation algorithm.

The corneal curvature is determined by measuring the distance between reflected light images as in conventional keratometry.

The anterior chamber depth is determined as the distance between the optical sections of the crystalline lens and the cornea produced by lateral slit illumination.

"White-To-White" is determined from the image of the iris.

Biometric eye measurements with the IOLMaster can be learned very quickly without time consuming training courses.

Extensive integrated safety features ensure both the patient and operator maximum safety when using the IOLMaster.

The software of the IOLMaster runs under Windows. The LC display functions both in patient eye alignment as well as a results and calculating interface. The instrument is controlled via the joystick and computer keyboard with integrated touchpad (mouse).

The integrated computer calculates the optimal intraocular lens on the basis of the measured values and international calculation formulas. Five formulas are integrated as standard.

Prior to use the data of the desired lenses must be entered into the database.

On the basis of postoperative refraction results, the lens constants that are entered in the calculation formulas may be individually optimized (personalized) for every user. In doing this, the data records of patients who had refractive surgery of the cornea should be ruled out.

Program control

General remarks

The control computer of the IOLMaster operates on the Windows operating system. For safety reasons, the operating system is not accessible to the user.

The known operating conventions of Windows analogously apply to the user interface of the IOLMaster software. This refers to the use of mouse/touchpad, the operation of buttons, the work with dialog boxes and menus, confirmation by double-click, etc.

Note:

The system does not support all key combinations of Windows. The special Windows keys that exist on some keyboards are ineffective.

The IOLMaster software uses only a few forced processes.
The user may optionally switch between the available modes.

In rare cases, it may happen that Windows error messages appear on the LC display. This might be the case, for instance, if the program run is affected (mostly by external disturbances).

Risks to persons are definitely precluded by multiply existing safety interlocks in hardware and software of the instrument.

If a warning message appears shut down the instrument and label it as being defective. Then, call Carl Zeiss Meditec Authorized Service.

The instrument may be controlled via:

- ☐ Buttons (using cursor, touchpad/mouse),
- ☐ keyboard or
- ☐ menus.

Some functions are alternatively operated by the pushbutton on the joystick.

Operation via touchpad and keyboard


Fig. 3 Operation via touchpad and keyboard

- ☐ Move the cursor by touching the touchpad with your finger and moving it as desired.
- ☐ Single and double clicks are possible by tapping a finger on the touchpad or with the left mouse button.
- ☐ To drag the cursor, hold the left mouse button depressed while moving the finger across the touchpad.
- ☐ The right mouse button is only functional for:
 - Resetting the zoom function (see p. 55)
 - Continuous fine positioning of the cursor while dragging (see p.59).

Single click	Selection of menu, textbox or entry. Operation of all buttons
Double click	OK, confirmation of actions.

- ☐ In addition to program control via mouse and touchpad, you may also activate particular menus by pressing individual keys or key combinations (see *menu overview on p. 17* and *overview of shortcut keys on p. 15*).

Operation via optional keyboard and mouse

You can also operate the instrument via standard PC keyboards and mice with PS/2 connector obtainable from computer shops in place of the keyboard with integrated touchpad provided. Connect these items to the appropriate connectors (see *Fig. 2., p. 10*).

Screen layout



Fig. 4 Screen layout (with software option A installed)

Overview of buttons and shortcut keys





Button	Key	Function	Remarks
	<N>	Activates data entry screen.	For inputting new patients as well as recalling patient data within the last 5 to 100 days Input of patient data is compulsory.
	<O>	Activates overview mode and surrounding field illumination.	For coarse alignment prior to measurement
	<A>	Activates axial length mode.	
	<K>	Activates corneal radius measurement mode.	
	<D>	Activates anterior chamber depth mode.	
	<W>	Activates WTW determination (optional)	WTW = <u>W</u> hite <u>T</u> o <u>W</u> hite
	<I>	Activates IOL calculation.	Calculation possible after measurement of one eye. After corneal refractive surgery, it is necessary to determine the corneal refraction.
	<P>	Prints results	
	<S>	Exports data.	Requires connection of an appropriate computer to the serial port. *
	<X>	Sends data to diskette.	Requires diskette to be inserted in disk drive.
	<E>	Exits IOLMaster software and Windows.	In case of problems: Turn off power switch and disconnect power cable!

* If connecting an external PC, the owner/operator must ensure to meet the safety requirements as per EN 60601-1-1 (medical electrical systems)!


Key	Function	Remarks
<Space>	Cyclic change of modes: ALM, KER, ACD, ALM, WTW	
Pushbutton on joystick	Program continuation Release of measurements	Effective in overview mode: Switches to ALM mode Effective in ALM, KER, ACD and WTW mode
	Deletes the selected ALM- or KER measurement from the list	Effective only in ALM-, KER- and WTW mode.
<Ctrl> + <D>	Effective in ALM mode: Prints the image of the selected graph. Effective in WTW mode: Prints the selected image of the eye.	ALM: One graph only. WTW: Right + left eye

Summary of result displays in Axial Length Mode

Display	Meaning	Explanations
	3rd axial length measurement (22.55 mm) selected.	The axial length graph of this measurement is displayed.
	Unreliable result $1.6 < \text{SNR} < 2.0$ SNR = Signal-to-Noise Ratio	"Borderline SNR" appears above axial length graph. The result should be examined by the user for validity.
	Erroneous result $\text{SNR} < 1.6$	"Error!" appears above axial length graph.
	Result has been manipulated.	* Remains displayed even if manipulation was undone!
	Signal-to-Noise-Ratio (SNR) and axial length next to axial length graph	Values of the peak under the measuring cursor.
	The measuring cursor is always placed above the signal peak with the highest absolute amplitude.	Confirmation of cursor correspondence to RPE signal peak is required.

Menu overview

The illustration below shows an overview of accessible menus and submenus for program operation via menu commands.

User Database

Since the instrument may be used by multiple surgeons at a group practice or clinic, surgeon specific records may be created. This is done via **User Database** on the **Options** menu.

- Click on **User Database** from the **Options** pull-down menu. This will bring up the dialog box for the entry of surgeon-specific data.

Fig. 5 Options menu

 A screenshot of a dialog box titled 'Please input password' with a close button (X) in the top right corner. The dialog contains several input fields: 'Name:' with a dropdown arrow, 'Password:' with a text box, 'New Password:' with a text box, and 'Comparison:' with a text box. There is a checkbox labeled 'Change Password' between the 'Password' and 'New Password' fields. On the right side of the dialog, there are two buttons: 'OK' and 'Cancel'.

Fig. 6 Dialog box for creating a surgeon profile

Note:

When the instrument is delivered, the User Database only contains the administrator without any password specifications.

Only the administrator is entitled to add or delete surgeons and edit their databases.

Warning:

Individual surgeons may edit their databases only if password protection has been set. If no password protection was set, the databases are accessible to all users!

The administrator may activate password protection after checking in the **Change Password** box.

- Type in the password in the **Password** and **Comparison** text boxes.
- Confirm your entries with **ok**.

- To create a new surgeon database the administrator opens his or her own database by selecting the **Administrator** option in the **Name** list box.

A dialog box appears where new surgeons may be added.



Fig. 7 User Manager for entry of new user and data backup

- Type in the name of the new surgeon.
- If several surgeons share the instrument it is recommended specifying a password each, which must be repeated in the **Comparison** text box.
- Use the **ADD** button if you have added a new surgeon or the **SET** button if you have changed the password of an existing surgeon.
- If you want to delete surgeon data from the database, click the **ERASE** button after having selected the name in the left window.
- Click **OK** to confirm your entries and log any new surgeons to the database.
- For the entry of lens data, refer to **Filling in the IOL database** (p. 37 f.).

Data backup

Backup (Creating a backup copy) With the Backup function, you can save the patient data used for the optimization of IOL constants and the IOL data of all operators and of the corresponding lenses used for the calculation to a diskette.

Independent of the set deletion date, the following patient data are saved: **Last Name, First Name, Date of Birth, ID Number, refraction data, visual acuity, remarks, examination date, mean value of axial length, mean value of corneal curvature** (the mean value used for IOL calculation), **mean value of anterior chamber depth** and **WTW data**.

The backup process also includes the tables used for IOL constant optimization (assignment of operator/lens/patient/eye/postop data). Additionally, for all operators, the IOL constants currently used for calculation will be saved.

Note:

In this way, all significant patient and IOL data and the data required for lens optimization can be saved.

Individual values of axial length, corneal curvature/refractive power, anterior chamber depth, WTW are not saved and may get lost, e.g. in the case of a hard disk fault.

You are therefore advised to regularly save these data to a diskette using the Export (X) function and thus minimize data loss in the event of a fault on the instrument.

Follow this procedure to create a backup copy:

- In the **User Manager activate Administrator**.
 - Click the **BACKUP** button to initiate the backup process.
 - Insert a formatted diskette in the disk drive.
 - Confirm with **OK**.
 - If the diskette should contain any files, you are asked to confirm overwriting of the files with **YES**. Answering the program query with **NO** aborts the backup process.
- While data is being copied to the diskette, a progress bar informs of the current state of the copy process.
- Finally, you will be informed that data backup was successful.

Note:

If the diskette should hold a previous backup copy of the database or other files, these will be deleted before the IOL database is copied.

Restore

Using the **Restore** function, you can retransfer saved data to the IOLMaster.


Follow this procedure to restore saved data:

- In the **User Manager** activate **Administrator**.
- Click the **RESTORE**.
- Insert the diskette with the latest backup copy; confirm with **OK**.
- A program query will ask you if all operators currently stored on the IOLMaster, the corresponding IOL data, and the patient data available for the optimization of the IOL constants can be overwritten.

On confirmation with **YES**, database data is being copied to the IOLMaster.

A progress bar informs of the current state of the copy process.

- Finally, the program will inform you if the restore action was a success.

Note:

After backed up data has been restored, the state of the user database in the **User Manager** corresponds to that of the time of backup.

Import

The Import function allows you to transfer IOL data (name and corresponding IOL constants) of a database saved to a diskette back to the IOLMaster. Imported data may be assigned to one or several operators. Import will only be successful, if the data on the diskette have the necessary file structure.

Follow this procedure to import IOL data:

- In the **User Manager activate Administrator.**
- Click the **IMPORT** button.
- Insert the diskette containing the database to be imported and confirm with **OK**.

Import


Fig. 8 Dialog box for the import of IOL data

- Click on "Lens constants <+>" to open the list of available lenses.
- Choose the desired lenses; select several lenses with <Ctrl> + cursor + click (selected lenses appear highlighted in blue).
- Choose the operator (one or several) with <Ctrl> + cursor + click (selected operators appear highlighted in blue); if not entered yet, enter new operators previously.
- Click on <>>> to copy the data.
A progress bar informs you of the state of the copy process.
Selected lens data will be added to the selected operators.
- Exit the dialog box with **EXIT**.

Setup submenu

- The **Options Setup** submenu provides the following settings:
 - Date/Time
 - Program Settings
 - Regional Settings
 - Keyboard
 - Add Printer
 - SW-Option
 - Update


Fig. 9 Setup-submenu

Date/Time

This option allows you to set Windows system date and time.


Fig. 10 Date Time

Program Settings/Program

This submenu provides the following options (defaults in bold types):

- Software language:
IOLMaster dialogs in German or **English** or other languages (change requires system restart).
- Display of visual acuity: Decimal or **Snellen**.
Entry of visual acuity in Patient Data dialog.
- Storage time of datasets: 5 ... **20** ... 100 days.
All figures between 5 and 100 are possible.
- Data Export:
Choose old, if the connected office management system only allows import of data of interface software versions 1.01 ... 2.02 (patient data, measured values).
Choose **new**, if the connected office management system can import all offered data according to interface software version 3.0 and higher.
- Keratometer display: Radius or **Corneal K's** for IOL calculation.
- Keratometer display: **one measure** oder **list of measures**.
By activating the "list of measurements" option, you may simultaneously display three measurements (joystick released 3x), consisting of five internal individual measurements each
You may take an optional number of individual measurements, but only the results of the last three measurements will be displayed.
- Keratometer display: – Cylinder or + **Cylinder**.


Fig. 11 Program Settings

- Keratometer refractive index:
Entry of equivalent refractive index for conversion of corneal radii into corneal K's.
Default: **n = 1.3375**.
Type in the refractive index implemented on your keratometer (refer to keratometer's User Manual).
- Printout of IOL calculation:
Choose if you want to have the calculated IOLs of **both eyes** printed on a single page or only one eye per page each.

Program Settings/Export (Option A)

On this submenu, you can choose the data to be saved to diskette:

- Patient ID:
Last Name, First Name, Date of Birth, Exam. Date, ID Number.
- Measured values:
Axial length (mean value, individual measured values),
Keratometer (current value; list of measured values),
Anterior chamber depth,
White-To-white (current value; list of measured values).

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Corresponding to Windows.

Keyboard

Change of the keyboard version corresponding to Windows.

Add Printer

Installation of a Windows printer.
Starts the Windows printer folder.

Software Option

Installation or uninstallation of a software version.

Update

Starts update installation of a new software version from a disk.

- Insert disk.
- Click on **Update**. The installation routine starts updating the software.
Follow the prompts appearing on the screen.
- Remove the installation disk from the drive.

Preparing measurements

Switching on

- Turn on the power switch. The IOLMaster starts by performing a self-test. Then, the data entry screen appears (Fig. 13).

Entry of patient data

Fig. 13 Patient data entry screen

The personal data of patients not listed yet in the database (New Patient) are to be entered via the keyboard. (Use "digits", "letters", ".", "-", and "_" only; other special characters are not accepted!).

- To move the text cursor to the next text box, use the <Tab> or <Enter> key or the mouse.

Note:

For safety reasons, it is absolutely necessary to type in surname, first name and date of birth. Data is saved as you type it (capital or small letters). The date of birth is accepted in the format mm/dd/yy and checked for plausibility. You may also enter the year as a four-digit figure (yyyy). With patients who are older than 100 years, this input format is absolutely necessary. The entry of visual acuity data is possible only in the data format set via **Options - Setup / Program Settings**.

- For the work with the database field, refer to p. 49 f.

- Conclude the entry by clicking the **<NEW>** button or striking the **<Enter>** key. This will automatically activate the "Overview" [OVW] mode. The fixation light and the illuminating LEDs are on.
- If requested, you may exit the application and Windows by a click on the Exit button.

Aligning the instrument to the patient

The two red ring marks (**13**, Fig. 1) on the side rails of the headrest serve for coarse vertical adjustment of the patient, (**3**, Fig. 1). The patient's eyes should be level with these marks.

- In overview mode align the instrument to the patient's eye using the joystick (**1**, Fig. 1). Vertical adjustment is by turning the control knob. Tell the patient to always look at the yellow fixation light in the middle.

Note:


In ALM mode, the fixation light is "red"; otherwise it is always "yellow".

- Adjust the instrument-to-patient distance until the 6 light spots (**1**, Fig. 14) appear focused.

Note:

The circle of lights should be approximately centered to the pupil and cross hairs (**2**, Fig. 14).

- 1** Circle of light spots for focusing
- 2** Cross hairs

Fig. 14 Video image of the eye with correctly aligned instrument

Axial length measurement [ALM]

- Activate the ALM mode using the
 - ALM button or the
 - <A> key or the
 - (pushbutton on the joystick).

Note:

With the ALM button and the <A> key you can start axial length measurements from all modes. The pushbutton on the joystick will work only in Overview mode.

To measure aphakic or pseudophakic eyes or eyes filled with silicone, select the corresponding mode from the **ALSettings** menu. This special AL mode is being displayed in the video image field and will be active until you reset it via the menu. The instrument will also be reset to the **"phakic"** mode by changing the side (other eye) or by measuring a new patient.

Note:

In the measurement of pseudophakic eyes and with specific intraocular lenses, two peaks may appear. The first higher peak is a side maximum of the IOL, while the second peak is produced by the retina. In this case, manual correction is necessary (see evaluation of ALM results).

On activation of ALM mode, the instrument automatically changes the magnification: a smaller section of the eye becomes visible with the reflection of the alignment light (**2**, Fig. 16) and a vertical line (**1**, Fig. 16).



Fig. 16 Video image of the eye with correctly aligned instrument

- Ask the patient to look at the red fixation light. On the display, a crosshair (**3**, Fig. 16) with a circle in the middle appears.
- Fine-align the instrument so that the reflection of the alignment light appears within the circle.

Key	Button
A; joystick push-button	


Fig. 15 ALSettings

- Start the measurement using the pushbutton on the joystick.

**Note:**

Before taking the measurement, tell the patient once more to fixate the red fixation light, which lights more intensively during the measurement. If the patient's visual acuity is poor, ask if she or he can see the fixation point! If not, erroneous measurements cannot be precluded.

In the case of poor visual acuity/high ametropia ($> \pm 6$ D), it is advisable to measure through the glasses. If taken conscientiously, such measurements will not produce erroneous results.

Measurements through contact lenses will lead to measurement errors and therefore should not be performed.

In the status bar, a blue progress bar appears. The corresponding display field beside the video image shows the measured axial length. The video image is overlaid with a graph similar to those known from ultrasonic measuring instruments. At the same time, the axial length and the signal-to-noise ratio (SNR) of the measuring signal are being displayed. The SNR is a measure of the quality of measurement. Measurements with an SNR between 1.6 and 1.9 appear with an exclamation mark after the measured value. At the same time, the message "Borderline SNR" appears on the display (Refer to the section **Post-run editing of axial length measurements**, p. 56 ff., for the assessment of the SNR).

With stronger lens opacities, it may be advisable to defocus the instrument. You may choose a reflection (2, Fig. 16) as large as the circle on the display. If even then measurements are impossible, refocus the instrument and shift the reflection to the bottom and/or top margin of the circle on the display by varying the vertical adjustment (turning joystick).

**Note:**

Defocusing and shifting the reflection within the circle is without any effect on the result, because interferometric axial length measurement is absolutely distance-independent.

If a valid measurement is achieved, you may assume that the displayed result is the correct axial length.

- By pressing the pushbutton on the joystick the next measurement of this eye is being started. It is possible to take up to 20 of such measurements per eye on one day.

**Caution:**

Avoid measurements of eyes with retinal detachment. In this case, erroneous measurements cannot be precluded.

After the second measurement the average of the axial length measurements will appear. This value will be recalculated and updated with each additional measurement. If the result of one measurement differs by more than 0.2 mm from the other(s), an "Evaluation" message will be displayed in place of the average. This indicates that the measurement results must be edited later (see *Post-run editing of axial length measurements, p.56 ff.*).

The last measured value is always highlighted by a blue background.

- You can immediately delete the current value if it is obviously the result of an erroneous measurement by pressing the **>Del>** key and confirming the action with the **<Enter>** key (see also *Post-run editing of axial length measurements, p.56*).

In the Mode field of the status bar beside "ALM" the number of measurements of the respective eye taken on this day is consecutively counted. If the count reaches 20, you cannot take further measurements of this eye on this day. It is also not possible to reset the counter. Deleted measured values (see above) do not affect the measurement count.


Fig. 17 Video image after axial length measurement

Measurement of corneal curvature [KER]

Key	Button
K; space bar	

- To activate the keratometer mode [KER], press:
 - the <space> bar (only after previous [ALM] mode or
 - the appropriate button or the
 - <K> key.

Tell the patient to fixate on the yellow light.

- Align the instrument so that the 6 peripheral measuring points are symmetrical with the circular crosshair and appear optimally focused. The central point usually is not focused. This point is not analyzed for keratometer measurement!

Note:

Be sure that all 6 peripheral points are visible and located in the field between the two auxiliary circles on the display. It is advisable to have the patient blink her/his eye shortly before the measurement to produce an adequate tear film. This will improve the reflectivity of the cornea. The measuring points should be circular or ellipsoid. If the measuring points are irregular (i. e. corneal scar) measurement is not possible.

Precise measurements are possible only if the 6 peripheral measuring points appear optimally focused on the display.

With dry eyes, possibly the administration of a tear replacement solution may help to obtain a correct measurement result.

- Start the measurement using pushbutton on the joystick.

Five measurements are taken within a period of 0.5 s and their average is displayed. The end of the measurements is indicated by a short acoustic signal.

The following is displayed: corneal curvature (in mm or diopters K) of the principal meridians with corresponding axes. If the cornea is spherical, only one radius or one refractive power value will be displayed.

On the display, a blue progress bar in the status bar indicates the process of computation.

The keratometer measurements may be repeated if desired.

If you take several measurements of corneal curvature in "one measure" mode, the previously measured values will be overwritten in the display.

To restore the last (just overwritten) results, press shortcut key <Ctrl> + <Z> (single-step "UNDO" function).

Note:

The "UNDO" function itself cannot be undone!


Fig. 18 Adjustment for corneal curvature measurement

Key
Ctrl + Z

Multiple keratometer measurements

If on the **Options – Setup / Program Settings** you choose the “one measure” option, only one measurement result will be displayed, which was obtained through five internal single measurements.

If on the **Options – Setup / Program Settings** you choose the “list of measures” option, three measurement results (joystick released 3x) may be displayed, obtained each through five internal single measurements. Though you may take an optional number of measurements, only the last three measurements will be displayed.


Fig. 19 Multiple keratometer measurement

- Same as before, you may delete the last measurement by pressing **<Ctrl> + <Z>**.
- To delete one of the three displayed measurements, select the measurement to be deleted and press the **** key. Confirm the delete action with “Yes”.

If the results of the last three measurements differ by > 0.5 D (average value of both radius) or the tolerance of $> 0.08 \dots 0.1$ mm (depending on n) has been exceeded, the display shows the message “Evaluation!”. In this case, check the tear film of the examined eye, let the patient blink, if necessary, and repeat the measurements. The erroneous measurements must be deleted as in the “Evaluation” state the obtained measured values will not be accepted for ACD measurement, IOL calculation and the database for optimization of constants.

- Highlight the result to be used for IOL calculation:
The last result will be highlighted in blue by default.


Fig. 20 Results with too large result differences

Anterior chamber depth measurement [ACD]

 **Note:**
The keratometer measurement must be performed before anterior chamber depth measurement!

Key	Button
D; space bar	

- To activate the ACD mode use the:
 - <**space**> bar (only from [KER] mode) or the
 - ACD button or the
 - <**D**> key.

The system will automatically turn on lateral slit illumination. This illumination subjectively appears to be very bright to patients. The patient should continue looking at the yellow fixation light.

- Fine-align the instrument so that
 - the image of the fixation point appears to be optimally sharp within the square on the display,
 - the image of the cornea is not disturbed by reflections,
 - the anterior crystalline lens is optimally visible.

As a rule, the image of the fixation point lies between the images of the cornea and the crystalline lens. It should be near (but not within) the optical section of the crystalline lens! The corneal image is unsharp due to system design.


Fig. 21 Alignment for anterior chamber depth measurement

 **Note:**
The alignment of the instrument, particularly in the case of small pupils, requires both operator practice as well as cooperativeness of the patient. The alignment procedure is easier on a dilated pupil (see also *"Tips on anterior chamber depth measurement"*, p. 71 ff.).

- Start the measurement using pushbutton on the joystick.

Note:

Before the measurement, tell the patient to look at the yellow fixation light and not into the lateral slit which is flickering during the measurement.

After engaging the measurement move the instrument to the side for patient comfort. The instrument calculates the ACD values, which takes approximately 15 seconds.

Explanation: Anterior chamber depth on the IOLMaster means the distance between the anterior vertex of the cornea and the anterior vertex of the eye lens. Hence, the displayed distance includes the thickness of the cornea.

Calculation of the anterior chamber depth requires the input of the corneal radius. If you performed a keratometer measurement before the measurement of the anterior chamber depth, the system will automatically use the measured radius for the calculation.

If the corneal curvature could not be measured with the IOLMaster, a window will appear requesting you to type in the radius (if the cornea is astigmatic, the values of both principal meridians).

- Enter a value between 4.0 and 13.0 mm (use decimal point). Proceed with **<OK>** or **<Enter>** key.

If you have chosen the display **"Refractive Index"**, enter a number between 26 and 80 please. When you enter the refractive power, make sure that on the IOLMaster the same keratometer refractive index is set as on the keratometer that was used for the measurement (--> p. 18).

In the display field beside the video image five ACD results are listed along with their average.

In the message bar, a blue status bar will appear. In the display field beside the video image, five ACD results are listed together with the calculated mean value.

The anterior chamber depth measurement may be repeated if desired.

If you take additional measurements of anterior chamber depth, the previously measured values will be overwritten.

To restore the last (just overwritten) results, press shortcut key **<Ctrl> + <Z>** (single-step "UNDO" function).

Note:

The "UNDO" function itself cannot be undone!

Key
Ctrl + Z

Determination of "White-To-White" [WTW] (optional)

Key	Button
space W	

- Activate the WTW mode as follows:
 <space> (only when in ACD mode) or
 Icon or
 <W> key.

Let the patient look at the yellow fixation light!
Align the instrument so that the six peripheral measuring points are symmetrical with the crosshair and the iris structures or the edge of the pupil appears optimally focused.

 **Note:**
Ask the patient, if she or he sees the fixation point. If the patient fails to fixate properly, the visual axis is not detected correctly, which may result in measurement errors.

- Start the measurement using pushbutton on the joystick.


Fig. 22 WTW Determination

The iris is imaged and its center determined. The WTW result presents the diameter of the circle segments visible on the screen as lying around the iris.

The image is evaluated by locating and determining the transition from the conjunctiva (white) to the cornea (colored iris).

Beside the WTW value, the deviation of the visual axis from the center of the iris is displayed. The respective values are given in millimeters with reference to a Cartesian coordinate system the origin of which is assumed to be in the found center of the iris. If the visual axis is above the iris center, the Y value will be positive; if it is below, values are negative. X values left of the center are negative; those right of the center are positive.

After you released the measurement on the joystick, one image of the eye is displayed. After having checked the correct detection of the iris, you must confirm the obtained results. Only then, the data can be printed and saved.

**Note:**

If the room lighting is inadequate, the software may fail to detect the iris or the fixation point sufficiently good.

It is advisable to keep extraneous light away from the instrument front panel and the examined eye.

You will obtain the best results with the examination room slightly darkened.

The WTW measurement may be repeated if desired.

Measurement of the other eye

The system registers automatically which eye is being measured (OD or OS).

All existing measured values of this patient are still in memory and may optionally be recalled.

- The measurement of the other eye is to be performed analogously to the previously measured eye.

**Note:**

After every changeover to the other side the overview mode is automatically activated for coarse alignment of the instrument.

Result printout

When the measurements are finished you may print the measured values and a measuring curve of that axial length measurement having the highest signal-to-noise ratio, a symbolic representation of the iris, the pupil and the WTW value. If you want to print the graph of a different axial length measurement, select it using the cursor (left mouse button) and press shortcut key <Ctrl> + <P>.

Note:
For operation of the printer, observe the user’s manual provided with the printer.


Fig. 23 Options menu

To set up an appropriate printer, click on the **Options** menu, **Setup**, **Add Printer**. This will bring up a dialog box for selecting and configuring the printer. The options available depend on printer hardware and software.

 **Note:**
Please, use printers recommended by Carl Zeiss Meditec only.

The printout of the measured values may be started from every measurement mode (ALM, KER, ACD, WTW). The printout will include all results obtained so far (also those of the other eye if already available).
It is advisable to start the printout only if all measured values of both eyes are available.

 **Note:**
While printing, do not take further measurements.

Key	Button
P	

- To start the printout:
 - click the printer button or
 - <P> key.
- To start the printout in WTW mode:
 - Icon or
 - <Ctrl> + <P> keys.

IOL calculation

If all measured values have been determined (depending on the calculation formula) the surgeon can use the IOLMaster to generate IOL options to match his/her surgical and postoperative requirements.

Filling in the IOL database

Before the system calculates IOL options, you must enter the desired lens types into the database.

- The user opens a database window for the entry of the specific lens data via the **Options - User Manager**. by entry/selection of her or his name (possibly, additional entry of password). (For adding a new user, refer to p. 16).


Fig. 25 Database window for entry of lens data

- In the **Name**, **A Const. Manufact** and **ACD Manufact** fields, type in the corresponding data of the manufacturer or the values specified in catalogues or on package slips.

Note:

If the ACD constant should not be available, you may click the **ADD** button after the entry of the A Const. Manufact. All parameters will automatically be calculated from the A Constant based on standard formulas.

- In the **A-Const. SRKII**, **A-Const. SRK/T**, **a0**, **a1**, **a2**, **pACD** and **SF** fields you must enter/edit your IOL constants optimized for different formulae of calculation for the optical biometry or the personally determined constants.

Note

For calculating with the IOL Master the suggestions for the power of the intraocular lenses to be implanted, do not use the IOL constants of the manufacturer, but constants optimized for the optical biometry only (see also p. 42 and p. 60).

- If you use lenses graded in 0.25 D intervals (in future), activate the **Power 1/4 D** radio button.
- Use the **ADD** button to save the data to the database.


Fig. 24 Input of password

- To delete the data of the lens type selected in the Lenses field, click the **ERASE** button.
- By clicking the **SET** button, existent data of a lens will be overwritten by any previously edited data.
- To enter the data of the next lens, overwrite the name of the lens. Exit the User Manager by a click on **OK**.

IOL calculation after refractive surgery of the cornea (optional)

Before calculating a proposal for an intraocular lens, the corneal refraction must be determined.

Taste	Icon
I	

- Initiate the calculation by:
Button (IOL calculation) or
<I> key.

Then:

- Select the tab "Prior refractive surgery".

Note:

This step of operation is necessary only with corneas pretreated by refractive surgery. With untreated corneas, IOL calculation is started instantly on selection of the biometric formula (see IOL calculation).

The corneal refraction is the quantitatively significant factor used in IOL calculation. Presently, it is impossible to exactly measure the corneal refraction that was subjected to corneal refractive surgery (e.g. by RK, PRK, LTK, Lasik or Lasek). That is why the corneal refraction must be determined in a different way. For this, two methods are available:

- ☐ Refractive History Method and
- ☐ Contact Lens Method.

For the **Refractive History Method**, the following values must be known:

- ☐ Preoperative corneal refraction (i.e. before the corneal refractive surgery,
- ☐ Preoperative refraction,
- ☐ Stable postoperative refraction and
- ☐ Corneal vertex distance.

As the change in refraction was achieved by the variation of the corneal refraction, the currently effective corneal refraction directly results from the difference between preoperative and postoperative refraction, corrected by the corneal vertex distance (vertex correction). The computational method is described in technical literature. If the corresponding data of the patient are available, the Refractive History Method delivers the most accurate results.

With the **Contact Lens Method** (contact lens overrefraction), one tries to determine the currently effective corneal refraction by two refraction measurements, once with and once without hard 'plane' contact lens.

The following parameters are needed:

- ☐ Refraction with contact lens,
- ☐ Refraction without contact lens,
- ☐ Refractive power of the (plane or almost plane) hard contact lens – refractive power of contact lens back surface and
- ☐ Corneal vertex distance.

In the ideal case, the refractive power of the contact lens back surface is equal to the unknown corneal refraction. From that point of view, several hard plane contact lenses with refractions of the back surface between 30 and 45 D should be available.

For the calculation of the corneal refraction, enter the necessary data of the patient in the data mask on the screen.

The screenshot shows the 'IOL Calculation' window with the 'SRK/T' tab selected. It displays data entry fields for both the right (OD) and left (OS) eyes. The 'Contact Lens Method' is selected for both eyes, and the 'Apply' checkbox is checked. The 'Clinical History Method' section shows pre and post-operative refraction data. The 'Contact Lens Method' section shows contact lens power, base curve, and refraction with and without the lens. The 'Corneal K's' section shows the calculated corneal curvature and radius.

Fig. 26 Correction of corneal refraction

For the calculation of the IOL, the Corneal K's selected by the examiner with <Apply> will be transferred to the IOL calculation table. On selection of the biometric formula, you can start IOL calculation.

IOL calculation

Key	Button
I	

- Start IOL calculation using the
 - appropriate button or the
 - <I> key.

The IOL calculation window appears. The measured values of both eyes will automatically be entered in the corresponding fields. The keratometer results are displayed in either Corneal K values (D) or Radius (mm) depending on your choice on the **Program Settings** menu under **Keratometer Display** (p. 18).


Fig. 27 IOL calculation and corneal refraction correction (after corneal refractive surgery)

- The IOL formulas **HofferQ**, **Holladay**, **Haigis**, **SRK/T®** and **SRK II®** are listed across the top.
Click on the appropriate tab to select the desired formula.
- By selecting his or her name from the Surgeon list box, the operator can access the operator-specific lens types saved to his or her data-base.
- The measured values may be edited if desired.
- Next, click on the eye whose IOL shall be calculated.
- Enter the desired target refraction. (No entry means plano).
- Select appropriate lenses from the available lens types.

- After you have entered the necessary data, click the **IOL CALCULATION** button to start IOL calculation for every selected lens type and every measured eye.

On the screen, only the data of the selected eye are being displayed. To view the calculated IOL data of the other eye, activate the radio button of the other Surgical Eye.

IOL Calculation

Fig. 28 IOL calculation

In the columns below the specified lenses, you will find the calculated refractive powers and target refractions for those lenses. The middle line appearing in bold types indicates that refraction of the corresponding IOL that comes closest to the desired target refraction.

Caution:

The IOL calculation is valid only, if the biometric measurement was correct, an appropriate IOL calculation formula was selected and the IOL constants were optimized before for the specific application.

- To print out the calculated data of the IOL to be implanted, click the **PRINT IOL CALCULATION DATA** button. The printout corresponds to the option selected on the **Program Settings** menu (page 18) showing the IOL data of both eyes on either a single page or on separate pages.
- The measured values and the previously selected IOL will be transferred to a correctly integrated PC (office management system) (up to software version 3.0, without Option A). Mark the IOL by a mouse click (mark: > ... <). You can mark only one IOL. The mark does not appear on the printout.
- With software version 3.0 plus Option A and the appropriate selection in the **Program Settings/Export** menu, all formulas, all selected lenses, all eyes and all proposed lenses will be transferred.
- Click on **OK** to finish IOL calculation.

Print IOL Calculation Data

Optimization of lens constants

The lens data available in the database, may be optimized by the following procedure.

- On the **Options** menu, activate **User Manager**. Select the respective operator and confirm your choice with **OK** to open the User Manager dialog box.
- Choose a lens. The input mask contains non-optimized constants calculated from "A Constant of Manufacturer" or constants optimized previously.

Fig. 29 Choice of lenses of a lens type

Optimize

- Click on the **OPTIMIZE** button.
In the dialog box of the selected lens, the lens constants will appear in the "Basis" column.


Fig. 30 Loaded lens type

- By a click on the **LOAD** button, the data records of all patients available for optimization will be loaded.

Fig. 31 Selectable data records

Special filter functions allow fast selection of patient data.
The right column shows the list of all patients available for optimization.

- Select the desired patient data record from this list.
- Choose the eye that shall be used for the optimization calculation. The fields below show the measured data of the IOLMaster.
- If you want the data of the other eye to be kept in the data table for the further optimization process, activate the **keep other side in table** check box (set check mark).
- Click on the << button to load the selected data record in the left-hand table. These data records are intended for IOL optimization.
- In this way, load at least 10 data records into the optimization table.
- If you want to retransfer a data record from the left-hand table to the right-hand table, select it and click the >> button.
- To delete a data record from either right-hand or left-hand table, select it and click on the button.
- If the left-hand table includes all data records you want to use, click on **OK**. The program will take you back to the optimization dialog box (Fig. 33).
- In following additional optimization processes, you may add further patient data record to the left-hand list.

Note:

The data contained in the database (right-hand and left-hand table) will not be deleted automatically. Thus, they are available also for later additional optimizations.

As hard-disk faults can never be completely precluded, it is advisable to regularly transfer the data to an office management system or print out the data.

Fig. 32 Selected data records


Fig. 33 Data records loaded for optimization

- Highlight a patient data record.
- In the **Impl. IOL (D)** text box, type in the power of the implanted IOL.
- In the **Post Op Ref** text box, type in the postoperative refraction.
- The entry of **Surgery Date** and **Post Op Date** is optional.
When entered, however, the data will be checked for plausibility.

Note:

Between **Surgery Date** and Post Op Date, there should be a period of at least 8 weeks. (This period, however, will not be checked!)

Note:

If a data record is highlighted in red, either no IOL and/or Post Op Ref data has been entered or a measured value (AL or KER) is missing for this data record!


Fig. 34 Data records prepared for optimization

- In this way, complete all selected data records. In the **Data Records** field, the number of data records containing IOL and Post Op Ref data (first figure) and the total number of loaded data records is specified.
- The text boxes below show the number of data records in the specified axial length ranges.

If you entered all IOL and Post Op Ref data, the requirements for the optimization calculation have been met.

Note:

If a data record is highlighted in yellow, it does not contain any ACD values yet:
a0 (Haigis formula) will not be optimized with such data records!

Note:

Data record appearing on a white background contain all data necessary for optimization.

- Start the optimization calculation by a click on the **OPTIMIZE** button. Depending on the number of data records to be processed, the computing process may take some seconds.

Next, the optimized lens constants will be displayed in the "new" column.

Note:

The displayed a0 value does not consider the data records highlighted in yellow!

The screenshot shows the IOLMaster software interface. On the left, a list of patient records is displayed, with the first record highlighted in yellow. The main window contains several input fields and buttons. The 'Last Name' field is set to 'Mustermann', 'First Name' to 'Wilhelm', 'Date of Birth' to '07/26/1916', 'Exam Date' to '08/14/2000', and 'ID Number' is empty. The 'Surgical Eye' is set to 'OS (left)'. The 'AL' field is set to '21.92'. The 'K1' field is set to '45.04' and 'K2' to '45.41'. The 'Opt. ACD' field is set to '3.1'. The 'Impl. IOL (D)' field is set to '20.5'. The 'Post Op Ref.' field is set to '0.75' and 'Post Op Date' to '09/18/2000'. The 'Surgery Date' is set to '08/15/2000'. The 'Data Records' field shows '46 / 46'. The 'Ave. AL' field is set to '23.53'. The 'Last Optimize' field is empty. The 'Basis' and 'New' columns are shown, with 'A Const.' values of '118' and '119.43', 'a0' values of '1.273' and '1.251', 'pACD' values of '4.96' and '5.808', and 'SF' values of '1.22' and '2.046'. The 'SRK II', 'SRK/T', 'Haigis', 'HofferQ', and 'Holladay' buttons are visible. The 'Load', 'New', 'Optimize', 'Erase', 'OK', and 'Cancel' buttons are at the bottom.

Fig. 35 Optimized lens constants

Fig. 36 Input mask for creation of a new data record

Caution:

In the fields for preoperative data, you may only enter data measured on the IOLMaster! When you enter the refractive power, make sure that on the IOLMaster the same keratometer refractive index is set as on the keratometer that was used for the measurement (--> p. 18). The entry of data measured on ultrasound devices will yield erroneous results!

- Complete the entries in the input mask.

Fig. 37 Filled-in data mask for a new data record

Note:

The entry of the measurement date is absolutely necessary!
Entry of **ACD data**, **Surgery Date** and **Post Op Date** is optional.

**Note:**

Between **Surgery Date** and **Post Op Date**, there should be a period of at least 8 weeks. (This period, however, will not be checked!)

- If you want to reject the entries made and return to the optimization calculation, click the **CANCEL** button.
- To confirm the new data record and add it to the list of data records to be used for optimization, click the **OK** button.
The new data record appears on the list of data records and is included in the **Data Records** field.

John, Q. Public, O.D. : lens # 1

new data record →

Mustermann, John, 12/05/1934, 07/05/2000

Mustermann, Antonia, 05/18/1921, 07/07/2000

Mustermann, Daniela, 06/25/1970, 08/17/2000

Mustermann, Eleonore, 12/11/1924, 08/04/2000

Mustermann, Elisabeth, 12/03/1915, 08/01/2000

Mustermann, Emil, 02/01/1934, 08/07/2000

Mustermann, Emma, 04/29/1922, 08/07/2000

Mustermann, Eva, 02/03/1920, 08/03/2000

Mustermann, Eva, 08/18/1920, 08/08/2000

Mustermann, Felix, 05/02/1909, 08/07/2000

Mustermann, Frieda, 01/11/1912, 08/03/2000

Mustermann, Harry, 07/31/1929, 08/11/2000

Mustermann, Hedwig, 11/22/1913, 08/04/2000

Mustermann, Hedwig, 11/29/1922, 08/14/2000

Mustermann, Helene, 02/02/1941, 08/01/2000

Mustermann, Helene, 04/15/1920, 08/04/2000

Mustermann, Helga, 02/03/1931, 08/16/2000

Mustermann, Hilde, 01/27/1921, 08/15/2000

Mustermann, Hilde, 02/15/1930, 08/04/2000

Mustermann, Hildegard, 05/05/1911, 08/10/2000

Mustermann, Ilse, 03/16/1927, 08/04/2000

Mustermann, Ingeborg, 09/25/1930, 08/14/2000

Mustermann, Johann, 07/13/1932, 08/02/2000

Mustermann, Käthe, 10/03/1914, 08/14/2000

Mustermann, Käthe, 11/18/1926, 08/01/2000

Mustermann, Klaus, 03/09/1926, 08/01/2000

Mustermann, Knut, 09/21/1970, 08/17/2000

Mustermann, Lillemor, 03/19/1901, 07/07/2000

Mustermann, Ludmilla, 08/05/1939, 07/04/2000

Mustermann, Manfred, 11/09/1929, 08/14/2000

Mustermann, Marga, 10/25/1929, 07/07/2000

Mustermann, Margarethe, 04/07/1915, 08/01/2000

Mustermann, Marianne, 09/20/1918, 08/08/2000

Last Name: Musterman Surgical Eye: OD (right)

First Name: John AL: 23.45

Date of Birth: 12/05/1934 K1: 43.25 @ 78

Exam. Date: 07/05/2000 K2: 42.75

ID Number: Opt. ACD: 3.6

Impl. IOL (D): 22.5

Post Op Ref.: -0.75 -0.5 x 90

Surgery Date: 07/08/2000

Post Op Date: 09/12/2000

Basis: << New: >>

A Const.: 119.13 << 119.13 SRK II

A Const.: 118.81 << 118.81 SRK/T

a0: 0.662 << 0.662 Haigis

pACD: 5.453 << 5.453 HofferQ

SF: 1.694 << 1.694 Holladay

Data Records: 47 / 47

14 .. 22mm: 6

22 .. 25mm: 35

25 .. 40mm: 6

Ave. AL: 23.52

Last Optimize: 09/18/2000

Load New

Optimize Erase OK Cancel

Fig. 38 New data record has been added

- Data records you do not want to use for a new optimization, but keep for a possible later reuse may be selected after a click on the **LOAD** button and retransferred to the result table by a click on the **>>** button.
If you click the **ERASE** button, however, the selected data records will be irrevocably removed.
- Click on the **OPTIMIZE** button to start a new optimization calculation.
- To reject the last optimization run, click on **CANCEL**. In this case, the optimized constants will not be saved to the lens data base although a new data record had been entered.
- Confirm the newly optimized lens constants by a click on the **<<** button right of the Basis field. In this case, all optimized constants will be accepted. If you want to accept a special constant only (e. g. a0), click on the **<<** button right of this constant.

- Click on **OK** to return to the **User Manager**.
The optimized lens constants will only be saved to the lens database and then be used for future IOL determinations if you confirmed them with **OK**!


Fig. 39 Optimized lens constants accepted

- Click on **OK** to return to the IOLMaster main module.

Note:

Data records with an IOL power of "0 D" will not be included in the optimization process.

In the optimization calculation, for every patient's data record those lens constants are calculated as they would have to be expected on grounds of the measured data and the results of surgery. Then, the average (sum of all lens constants divided by the number of patients) and the standard deviation are calculated.

The program will delete those lens constants that show a difference to the average larger than double the standard deviation. The thus resulting average will be displayed as optimized constant.

To obtain optimum constants, you should exclude patients with preoperative, intraoperative or postoperative complications that might affect the refraction state.

New patient

To add a new patient to the database, use the

- appropriate button or the
- <N> key.

Key	Button
N	

The results of the previous patient displayed on the right and/or left are cleared from the display. They remain however stored in the internal database for 5 to 100 days (—> **Program Settings**, p. 23). On the screen, the dialog box for patient data entry reappears.

Note:

The patient data and measurement results will not be deleted from the data tables for the optimization of lens constants (Fig. 31 and Fig. 33).

- After you have typed in and confirmed the data of the next patient, the program will switch again to the overview field with surrounding field illumination to assist in patient alignment.

Note:

The above order of measurements is only an example. You may also run the measurements described above in a different order. The only requirement is that the keratometer measurement precede the anterior chamber depth measurement.

Using the patient database

text box
for search word


Fig. 40 Patient data entry screen

The IOLMaster keeps an internal patient file for a limited time. The data is saved there and recallable (for reviewing, editing, printing and exporting data via serial port, export to diskette) for a period of 5 to 100 days (—> *Program Settings*, p. 18).

Note:


The IOLMaster is not designed for archiving patient data and measurement results.

The database field is structured like Windows Explorer. A "+" sign at the branch indicates that the database already holds measurement results for this patient.

- If you click on the "+" sign the data of the previous measurement(s) becomes visible.
If you click on the "-" sign the branch is collapsing again.

The data records of the file are sorted alphabetically by last name. To quickly access a data record, you can use the Search textbox. If you place the cursor in this box and type in the desired last name, all relevant data records will be listed.


If a patient visits the office again before the selected interval for deletion has reached, his/her data can be instantly transferred to the input area by clicking the patient's name.

- If you want to take a new measurement, click the **<NEW>** button or the shortcut key **<Alt> + <N>**.

Recalling previous measurement results

The system permits data records of previous measurements to be reviewed.

- Click on the "+" sign in front of the patient's name.
- Use the cursor to select the date of interest.
- To recall the data record, click the **OPEN** button or the shortcut key **<Alt> + <O>** or double-click on the date.

You can view and edit the data record.

- Automatic Right/Left detection is deactivated.
To choose a side, press the **R** or **L** key, or click with the mouse onto the desired side.

Erasing a patient

To erase a patient from the list, choose one of the following options:

- Select the patient (by a click on the name) and press the key or
- on the **Patient** menu, activate the **Erase** function.
- Confirm the erase action with **YES**.


Fig. 41 "Patient" menu

Renaming a patient

If you want to edit the Last Name, First Name, Date of Birth or the ID Number of a patient, follow this procedure:

- Select the patient (by a click on the name). Then, on the **Patient** menu, activate the **Rename** function or press <Ctrl> + <U>.



Fig. 42 Renaming a patient

- In the appearing dialog box with the patient data, edit the desired entries.
- On confirmation by a click on the **RENAME** button, the data of the patient will be changed in all relevant files.

Data export (optional)

Data may be exported to connected office management systems, personal computers or diskettes.

Data export to office management systems can only be implemented by the supplier of such systems. Please contact the supplier. Before the export of data, select the export mode appropriate for the connected office management system via Options → Setup → Program Settings.

For data export to a connected personal computer with a Windows operating system (Windows 95/98/2000/NT/ME), Carl Zeiss Meditec supplies the corresponding accessories. These include a serial cable (null modem, female/female connector) and software (on CD-ROM) to be installed on the PC. Data is stored to a database on the PC. From there, data may be exported to other file formats. The graphs of axial length measurements are made available in JPEG format.

Data export includes patient data, measured values and calculated IOL's. To activate data export to an office management system or a PC, click the

- corresponding button or
- press the <S> key.

Data is transmitted through the serial cable to the connected PC.

Note regarding data export to a PC:

- The PC must have been switched on and the software for data receipt started.
- On the screen of the IOLMaster, an export progress bar is visible.
- Data can be archived on the PC or processed appropriately.

Patient data may also be exported from the office management system to the IOLMaster.

For data export to a diskette:

- Click the corresponding button or:
- Press the <X> key.

Data is made available in a text file (IOLMexport.csv) for archiving and data analysis.

Key	Button
S	

Key	Button
X	

Switching off the IOLMaster

Note:

The procedure described below does not apply if there is a problem or if the instrument does not respond to your entries!
If this occurs, immediately switch off the instrument and disconnect the power cable. Label the instrument as being defective and call Carl Zeiss Meditec Technical Service.

- After you have performed all measurements and calculations, exit the program using the
 - appropriate button or the
 - **<E>** key.
- Confirm the operation with **<OK>** or **<Enter>**.
The data of the current (last) patient will be saved automatically.
- When Windows has shut down (wait for corresponding message!), you may switch off the instrument with the power switch (**1**, Fig. 2).

Key	Button
E	


Note:


Failure to switch off the device as described (emergency) may result in data loss and software damage.
To restart the device, observe the instructions given in Section *What to do in case of problems or system errors*, p. 81.

Example signal curves of axial length measurements


Valid signal curves

- Very good signals (signal-to-noise ratio > 10)
- Several secondary maxima visible (system-specific)
- Clear media, correctly fixating patient
- Weak ametropia

- Clear signal (SNR > 2.0)
- Secondary maxima visible
- Relatively clear media


- Signal in "Borderline SNR" range of 1.6 ... 2.0
- Steep rise of measuring signal
- On the display, results of such measurements are marked by an exclamation mark. Additionally, the message "Borderline SNR" appears.

This measurement result may be used after verification and comparison with other data of this measurement series.

Non valid signal curves

- Low signal (signal-to-noise ratio < 1.6)
- Error message is displayed.
- The measuring signal cannot be clearly distinguished from the noise.

Possible reasons:

- unsteady (non fixating) patient
- strong ametropia
- dense medial opacity along the visual axis


Repeat the measurement!

Ask the patient to fixate steadily.

Zooming the graph display

The system allows zooming the graphs in 4 steps for improving the interpretation of signal curves:

- Move the cursor on the length axis (X-axis) to the position that shall be in the center of the zoomed in image and press the left mouse button. You may repeat this procedure four times.
- To return to the original view (zooming out), place the cursor at any position on the length axis and press the right mouse button.

Note:
In zoomed views, the axial length scale is not visible.

Post-run editing of axial length measurements

The results of axial length measurements must be interpreted on the basis of the signal-to-noise ratio and the appearance of the graphs (cf. *"Example signal curves of axial length measurements"*, p.53). The illustrations below do not contain the video image.


Fig. 43 Presentation of the graph of the fifth axial length measurement (23.00 mm; SNR: 6.6) on the right eye (shown without video image)

SNR categories

While the system is internally calculating the axial length from the interference signal, it automatically analyzes the SNR.

SNR > 2.0 --> The measured value is valid.

SNR in the range 1.6 ... 2.0 --> The measured value is not reliable
(Fig. 44)

Fig. 44 Non reliable value

The signal-to-noise ratio may be low for the following reasons:

- Dense medial opacity along the visual axis,
- restless patients,
- alignment of instrument to patient eye is not optimal,
- very high ametropia (> 6 D),
- corneal scars,
- retinal pathology.

Note:

In this case, "Borderline SNR" does not mean a wrong result, but shall only remind you to verify this measurement!

Accurate signal peaks can be determined by comparison with other values of the measurement series of this eye (and of the other eye, if necessary). See also the sections "Example signal curves of axial length measurements", p. 53, and "Shifting the measuring cursor" below.

SNR < 1.6 (The measured value is unusable)

It is marked on the display as an erroneous result (Fig. 45).

This implies that the true measuring signal does not clearly stand out from noise. As a rule, the results of such measurements are not usable and should be deleted. You may delete them from the list by using the key once they are highlighted.


Fig. 45 Erroneous result

Note:

Erroneous results (display: Error) are not considered in mean value calculation!

The measured value may be transferred to the list of measured values by a click on the measuring cursor (white dot).

See to it that the measured values are consistent.

Shifting the measuring cursor

The measuring cursor (white point) is automatically placed on the center of the signal peak with the highest absolute amplitude. The corresponding axial length value is displayed beside the graph and in the display field. The SNR is calculated and displayed for this signal peak. The measuring cursor is placed in the center between the regions corresponding to half the maximal amplitude. If the signal curve is symmetrical (Gaussian curve), the cursor is positioned exactly above the maximum of the signal. For interpretation of axial length measurements (p. 60) it may be necessary to shift the measuring cursor to verify axial length agreement between values in a measurement series.

There are two ways to move the measuring cursor to another peak. In both cases, it is advisable to do this with the signal curve being zoomed:

1. Automatic positioning on a "remote" peak:

- Place the mouse (arrow) cursor in the measuring cursor (white point). Holding the left mouse button depressed move the measuring cursor over the desired peak. A vertical blue line will appear below the measuring cursor for easier orientation. This line will be moved along with the cursor.
- When you release the mouse button, the measuring cursor will automatically "snap in" over the desired peak.

Fig. 46 Value denoting manipulation

Then, the display will show the corresponding axial length value and the corresponding SNR, which will always be smaller than the maximal SNR found automatically. This recalculated measured value will appear in the display field denoted by an asterisk (*) (Fig. 46).
Example:


Fig. 47 Moving the measuring cursor to a different peak (signal curve zoomed in 3-times)

Explanation:
This manipulation will work only if the measuring cursor is moved fully across the (local) maximum of the desired target peak. This procedure is necessary for the search algorithm to reliably find the desired peak without returning and snapping in to the original (higher) peak. Closely adjacent peaks (double peaks) cannot be separated by this automatic method unless the curve adjoining them drops down below a value which is less than half the amplitude of their maxima.

 **Note:**
While dragging the measuring cursor, the original axial length value and SNR remain displayed beside the signal curve. The new axial length value and the corresponding SNR will be calculated and displayed only when you release the mouse button.

2. Fine-shifting the measuring cursor

Note:

This manipulation always requires previous zooming in!

- Proceed as described above under Item 1, but use the right mouse button for dragging the measuring cursor. This way the automatic peak detection is deactivated and the white point may be positioned at any place over the measuring curve.
- When you release the mouse button, the current axial length and the new SNR will be calculated and displayed.

This kind of manipulation is advisable with closely adjacent double or triple peaks.

Also in this case the recalculated axial length is denoted in the display field by an asterisk (*).

Explanation:

Even if you undo the manipulations with the measuring cursor (by moving it back to the automatically found maximal peak) and the measured value agrees with the original one, the asterisk after the measured value will remain. This is to indicate that the curve has been manipulated!

The described manipulations of the measuring cursor may be performed both in axial length measurement mode (after the current individual measurement) and in post measurement editing mode.

Interpretation of axial length measurements

As a rule, an interference signal is produced if the measuring light is reflected by the tear film and the retinal pigmented epithelium of the eye. This signal is utilized for axial length measurements.

Note:

Ultrasonic biometrical instruments measure the axial length as that distance between the cornea and the inner limiting membrane because the sound waves are reflected at this membrane. To ensure that the measured values obtained with the IOLMaster are compatible with those obtained through acoustic axial length measurement, the system automatically adjusts for the distance difference between the inner limiting membrane and the pigmented epithelium. The displayed axial length values are thus directly comparable to those obtained by immersion ultrasound!

It is absolutely necessary to re-personalize the "lens constants" for use with the IOLMaster prior to applying its calculation values to determine lens power and surgical technique. Refer to the specialist literature and the publications of the authors of the IOL formulas regarding the personalization of constants (Surgeon Factors).

Updated information is available on the Internet at:

http://www.meditec.zeiss.com/iol_master and/or

<http://www.augenklinik.uni-wuerzburg.de/eulib/>

With an optimally aligned instrument, a good SNR and weak ametropia (approx. ≤ 6 D), the secondary maxima located symmetrically on each side of the measuring peak will be seen (Fig. 48).

These maxima are an artefact of the light source.

The distance from each of the secondary maxima to the peak signal is approximately 0.8 mm. The secondary maxima are likewise always visible in measurements of the provided test eye.


Fig. 48 Undisturbed measuring signal with secondary maxima


The measuring system of the IOLMaster is capable of resolving fine structures on the fundus of the eye.

Depending on the anatomical conditions of the measured eye it may happen that the measuring beam produces interferences when being reflected at the inner limiting membrane and/or the choroid.

This is indicated by

- broader signal peaks of the measuring curve,
- variations of axial length data of approx. 0.15 - 0.35 mm in one measurement series and
- display of "Evaluation" in place of mean value.

Example:

Such measuring curves or measurement series absolutely require verification. This can be done both between the individual measurements (in ALM mode) and in post measurement editing (without the patient sitting in front). Interpretation or post measurement editing should always be done with the help of the zoom function.

Note:

The resolution of fine retinal structures is clearly distinguishable from the previously mentioned secondary maxima. The secondary maxima (2) are much further distant from the retinal pigment epithelium signal peak and they are symmetrical to it.

The distance between the maximal peak and internal limiting membrane or choroid is 0.15 - 0.35 mm (whereas the secondary maxima are about 0.8 mm displaced from the maximal peak!).

Signals from the inner limiting membrane (ILM)

Measuring light is relatively often reflected at the inner limiting membrane thus also producing an interference signal. The corresponding ILM signal peak appears to the left of the measuring peak for the RPE by a distance of 0.15 - 0.35 mm. If the measuring cursor is placed at the ILM it will produce a shorter inaccurate axial length value by 0.15 to 0.35 mm. This is the most likely reason that no average axial length value is generated for a series. In the zoomed view of the measuring curve, both peaks may be clearly distinguished from each other (Fig. 49).


Fig. 49 Double peak produced at inner limiting membrane (left) and RPE (right) (triple zoom)

Usually, the signal from the amplitude of the peak produced by the inner limiting membrane is smaller than that of the reflectance from the pigmented epithelium. In such a case the automatic algorithm finds the correct axial length.

Note:

Never move the measuring cursor manually to the (left) peak produced by the inner limiting membrane!

In rare cases it may happen that the amplitude of the signal from the inner limiting membrane is higher than that of the reflected light from the pigmented epithelium. In this case automatic peak detection will place the cursor incorrectly at the peak from the ILM.


Fig. 50 Signal curve with higher signal from inner limiting membrane (left peak) (double zoom)

In measurement series such individual measurements stand out by deviations in the range of approx. 0.15 - 0.35 mm towards shorter axial lengths. You may correct the measured value by moving the measuring cursor to the right to the smaller peak (that produced by the pigmented epithelium). This manipulation is permissible only with the other signal curves of this measurement series!

Signals from the choroid

Triple peaks

In rare cases it may happen that the measuring light is also reflected from the vessels of the choroid.


Fig. 51 Rare triple peak (triple zoom)

The measuring peak produced by the choroid appears shifted towards longer axial lengths by approximately 0.15 - 0.25 mm from the peak of pigmented epithelium.

In the example illustrated above, the signal from the RPE (middle peak) has the highest amplitude. The automatic peak detection system has correctly placed the cursor here to detect the correct axial length. No further movement of the cursor is needed. This type of rare triple peak clearly differs from the secondary maxima produced through the light source by the distance from the RPE reflected peak (Fig. 51).

In rare cases, depending on the anatomical conditions of the measured eye, it may occur that the signal produced by the pigmented epithelium does not have the highest amplitude.


Fig. 52 Triple peak (double zoom)

In Fig. 52 the automatic peak detection system will find an axial length value that is too short by approximately 0.15 - 0.35 mm. After comparing all measured values and curves for this eye, the measuring cursor must be moved manually to the middle (smaller RPE) peak produced by the RPE. This measured value may be corrected (with corresponding asterisk) or deleted.

Double peaks

In very rare cases it may occur that signals are produced by both the pigmented epithelium and the choroid.


Fig. 53 Double peak produced by pigmented epithelium and choroid (double zoom)

In Fig. 53, the automatic peak detection system has placed the measuring cursor in the correct position at the RPE and the measuring cursor must not be shifted manually.

Note:

Such a curve may only be evaluated correctly by viewing all measuring curves of this eye. Such a curve must be clearly distinguished from double peaks produced by the inner limiting membrane and the RPE (see Fig. 50)! It may be necessary to delete suspicious curves and take additional measurements (20 are maximally possible per eye on one day).

How to adjust the measuring marks

Ask the patient to look at the yellow fixation light. If the patient cannot see the fixation light instruct him or her to look straight ahead into the instrument.

Note:

The peripheral infrared measuring marks are invisible to the patient. (However in a darkened room an attentive observer may perceive the measuring marks as faint red dots when looking into the projectors of the keratometer.)

When aligning the instrument, take care that **all the six peripheral measuring marks** are visible and that they are lying as concentric as possible to the two auxiliary circles on the display. By varying the distance between patient and instrument, **optimally focus** the images of the six peripheral measuring marks on the display. The images of the measuring marks should be circular or ellipsoid.

To improve the reflectivity of the cornea it is advisable to ask the patient to close and open the eyes several times. In this way, the tear film is replenished and imaging of the measuring marks improved.

As a reminder, an appropriate message will appear below the video image when you activate the keratometer mode.


Fig. 54 Optimally aligned instrument (without crosshairs and auxiliary circles; the central fixation point is distinctly fainter than the six peripheral measuring points)

Explanation:

Depending on the reflectivity of the cornea the display monitor image of the central fixation point may be barely visible. This is insignificant for the calculation of the corneal curvature. Most important is sharp focusing of the peripheral marks.

Measurement errors

The "Error" message may have two basic causes:

- ❑ The measured values of the internal individual measurements vary by more than 0.05 mm (very rare or, if instrument was defocused).
- ❑ The measuring marks are either not recognized or recognized only partly as such. (The measuring marks not recognized will be identified on the screen after the measurement).

The possible reasons for this are described below:

Maladjustment

Defocused measuring marks


Fig. 55 Image with defocused instrument

Cause

The images of the measuring marks are too large due to the instrument being defocused. The system cannot calculate a measured value. The display field shows "Error".

Solution

After correcting the focus adjustment to minimize the peripheral mark size the measurement can be retaken (old value will be automatically overwritten). Sometimes, with exactly adjusted focus, small circles (like haloes) may be visible around the six peripheral measuring points. In this case, focussing is optimal.

Missing measuring marks


Fig. 56 The upper two measuring marks are blocked by the eyelid

The display field shows the "Error" message. This error may also occur, if the patient blinks his or her eye during the measurement (0.5 s). This may happen particularly with noncompliant or anxious patients.

Ask the patient to open his or her eyes wide and repeat the measurement.

If this does not resolve the problem carefully pull up the upper eyelid, like as known by tonometry, without indenting the globe.

Pushing against the globe will produce an erroneous measurement.

Cause

Solution

Other findings

Pseudophakic eyes


Fig. 57 Pseudophakic eyes

Cause

In the measurement of pseudophakic eyes, it may happen that images produced at the front side of the intraocular lens are visible beside the reflections from the cornea.

The reflections from the IOL are fainter and unsharp.

Solution

If no measurement is achieved after focussing the peripheral measuring marks from the cornea move the instrument approximately 1 mm away from the patient's eye (defocusing) and release the measurement. The images produced at the cornea will be slightly larger, while the artifacts of the IOL become fainter such that the evaluation process possibly does not identify the IOL reflections as measuring points. If this procedure does not succeed you cannot measure the corneal curvature.

Dry eye


Fig. 58 Multiple reflections produced by a dry eye


Fig. 59 Light trail (bottom) due to a dry eye (on top additional disturbance by an eyelash)

If the tear film is quickly broken up the reflectivity of the cornea will be greatly reduced at these points. If a measuring mark is projected to such an area, the otherwise circular or ellipsoid image of the measuring mark will become irregular. You will see irregular marks and/or multiple reflections. In this case, a precise measurement of the corneal curvature will not be possible. The results will fluctuate or the "Error" message will be displayed.

Cause

Have the patient close and then open eyes several times to replenish the tear film on the cornea and take the measurement immediately afterwards.

Tear film supplements may also be used to aid in measurement.

Solution

CauseIrregularities of the corneal surface (scars)

Scars and local irregularities on the corneal surface impair the imaging quality of the measuring marks. Depending on the expansion and location of these irregularities erroneous measurements may arise.


Fig. 60 Local corneal scar impairing measuring points at 3 and 5 o'clock

Solution

Try to position the measuring mark beside, above or below the scar by slightly displacing the instrument relative to the eye. Then, release the measurement. In such cases, it is advisable to repeat the measurement several times. Depending on the degree of irregularity fluctuating or erroneous results may occur.

Fig. 61 Condition after keratoplasty

In Fig. 61, it is impossible to take a keratometer measurement using the IOLMaster.

How to adjust the instrument

Ask the patient to look at the yellow fixation light at all times. If the patient cannot see the fixation light, let him or her look straight ahead into the instrument. When you turn on the anterior chamber depth mode, the system automatically activates lateral slit illumination. The illumination is always coming from a temporal direction to the eye. The slit illumination will subjectively appear bright to the patient. The measured values of the light load, however, (see *Technical Data*, p. 88) is smaller by some orders of magnitude compared to slit lamp examinations.

When you release the measurement, the lateral slit illumination starts flickering.

(Note: Looking into the slit projector is under no circumstances dangerous but it will lead to erroneous anterior chamber depth values.)


Fig. 62 Optimally adjusted optical section for anterior-chamber depth measurement

On the display, an image similar to that of a slit lamp is visible. It represents an optical section through the anterior segment of the eye. Align the instrument to the patient's eye by lateral adjustment using the joystick until:

- ☐ The image of the fixation point appears optimally focused in the square on the display,
- ☐ the image of the cornea (outer or temporal of the two optical sections) is relatively free of reflections (although unsharp is ok) and
- ☐ the image of the anterior crystalline lens is visible in the pupil.

Note:

The image of the fixation point should be near (but not in!) the image of the lens.

Tips for anterior chamber depth measurement

If the instrument has been properly aligned, the images of the fixation point and the front surface of the crystalline lens will be sharp at the same time, as they are approximately in the same plane.

As a rule, the image of the fixation point lies between the image of the anterior lens and that of the cornea if the instrument is optimally aligned.

**Note:**

The image of the fixation point must not be placed in the image of the lens or cornea!


Fig. 63 Optimally adjusted optical section (lens with cataract)

Fig. 62 and Fig. 63 show optical sections of right eyes.

The patterns left of the corneal image are direct reflections of the luminous light exit aperture of the lateral slit projector. These reflections are not needed for the calculation of the anterior chamber depth. They must **not** affect the image of the cornea (see below).

At the left margin of the picture, additionally reflections are visible produced by the ambience (here, a window) being in front of the patient. Depending on the lighting conditions in the examination room it may also be possible that you see the front side of the IOLMaster produced by the cornea. These artifacts do not disturb the measurement of anterior chamber depth unless the significant image details (images of cornea and crystalline lens) as well as the image of the fixation point are significantly weaker than some kind of extraneous light. This may be alleviated by slightly darkening the examination room.

If the above requirements are not met the results will be either erroneous or not possible. Because of the complexity of the images measured it might be possible that erroneous measurements are not recognized as such under certain circumstances.

The IOLMaster must be adjusted very carefully for anterior chamber depth measurements.

The measurement of the anterior chamber depth on eyes with very small pupils (e.g. with glaucoma) is particularly problematic and needs some practice.

The anterior chamber depth of the human eye also depends on the accommodative state of the eye. This cannot be assessed from an optical section of the anterior segment.

**Note:**

It is advisable to measure accommodating patients under cycloplegia.

Measurement errors

Like in keratometer measurements, there may be two basic reasons giving rise to the display of the "Error" message:

- ☐ The results of the five internal individual measurements vary by more than 0.15 mm (very rare) or
- ☐ the produced images (optical sections) do not contain relevant structures (normally, the edge of the crystalline lens) or disturbances are preventing their detection.

Eyes that cannot be measured:

Aphakic eyes cannot be measured, as the calculation of the anterior chamber depth is based on the evaluation (distance measurement) of the image of the optical section.


Fig. 64 Optical section of an aphakic eye; no lens visible

Tips for anterior chamber depth measurement

The same refers to pseudophakic eyes if the intraocular lens does not sufficiently scatter the light from the slit.


Fig. 65 Optical section of a pseudophakic eye. On the right, the unsharp reflection of the rear face of the IOL is visible.

Under certain circumstances, a strong reflection may be visible on pseudophakic eyes, which is produced by the rear surface of the lens.

When you release an anterior-chamber depth measurement in such a case, the evaluation software will interpret the bright reflection as the front surface of the crystalline lens. The instrument will then display a measured value in the range of about 5 to 6 mm. **This value is not the exact anterior chamber depth!**

Maladjustment

Defocused instrument


Fig. 66 Image with strongly defocused instrument

If the instrument is not optimally focused, the image of the fixation point becomes larger and fainter. At the same time, the images of the front edge of the lens and/or the cornea may become so faint that they cannot be recognized as such. In such a case, the system displays an "Error" message and an explanatory text indicates which image details are either missing or could not be recognized correctly.

Cause

Improve the focus adjustment of the instrument and repeat the measurement. The fixation point must be optimally focused.

Solution

Note:

As a rule, slight defocusing of the instrument does not affect the anterior chamber depth measurement significantly.

Missing lens image with phakic eyes

Particularly in the case of eyes with small pupils it may happen that no light is reflected back into the viewing optics of the instrument. A slight lateral misalignment may make the lens invisible. This problem may also appear with restless or poorly fixating patients.

Tips for anterior chamber depth measurement


Fig. 67 Slit image on the iris (invisible lens)

Cause

In such a case, the slit image on the iris is (almost) continuously visible. The automatic evaluation software does not recognize this kind of maladjustment. The system will display values that are too short. These values do **not** correspond to the actual anterior chamber depth, but represent the distance between the anterior cornea and the iris. Make sure that the lens section is placed on the front section of the crystalline lens.

Solution

Adjust the instrument laterally until the anterior lens becomes visible. If necessary, ask the patient to look steadily at the fixation light. Then, repeat the measurement.

Note:

It will suffice if a relatively small section of the lens is visible. The picture below shows an alignment which permits accurate measurement.


Fig. 68 Minimally visible anterior lens

This image is sufficient for the calculation of the anterior chamber depth. (In this photo, the front side of the IOLMaster is visible as a non-disturbing artifact).

Image of fixation point in image of crystalline lens


Fig. 69 Fixation point lying in lens image

If the image is laterally misaligned, it may happen that the image of the fixation point is lying within the lens image.

Cause

Position the instrument so that the fixation point lies between the images of crystalline lens and cornea. Then, repeat the measurement.

Solution

Reflections in the corneal image


Fig. 70 Reflection in corneal image due to lateral maladjustment

The lateral adjustment of the instrument is not correct. In the corneal image slit illumination reflections can be seen and the fixation point lies within the lens image.

Cause

Adjust the instrument laterally until the corneal image is undisturbed. As a rule, the fixation point will then be between the image of the anterior lens and that of the cornea. Repeat the measurement.

Solution

Pathological findings

Dry eye



Fig. 71 Optical sections of dry eyes

Cause

A locally disturbed tear film changes the scattering properties of the cornea considerably. For this reason, the optical section of the cornea may become irregular.

Solution

Let the patient close and open their eyes several times to replenish the tear film on the cornea. Release the measurement immediately afterwards or use a tear supplement to prevent quick tear breakup time.

Irregularities of the corneal surface (scars)


Fig. 72 Condition after keratoplasty (same eye as shown in Fig. 61, p. 70, section *Tips for keratometer measurements*)

Scars and local irregularities of the anterior cornea impair the image quality of the optical section of the cornea. Depending on the expansion and the degree of these irregularities, this may lead to measurement errors.

Cause

In such a case, the fixation point is imaged as a "cloud" and it is impossible to improve the adjustment. If nevertheless apparently plausible anterior chamber depth results are displayed, they can only be regarded as reference points.

Solution

To obtain reliable data, all the other well known facts and findings of this eye should be included in the evaluation.

How to adjust the instrument

Ask the patient to look relaxed at the yellow fixation light.

Focus on the iris, not on the peripheral points. Adequate room lighting will facilitate detection of iris structures. Avoid direct exposure of eye and instrument front panel to extraneous light.

Particularly take care that the visible right and left edge of the iris is not disturbed by reflections from lamps and windows.

If the iris structure should not be discernible, focus on the edge of either iris or pupil.

Larger degrees of defocussing will result in incorrect data.

After the image has been taken, the operator should check if the software correctly detected the edge of the iris. If the circle segments drawn in the image should not define the iris correctly, discard the result. Click on **OK** to confirm the results and save the data.

What to do in case of problems or system errors?

To restart the instrument after a system error or improper exiting of the program.

- Turn power switch to I.
 - The Scandisk utility will automatically run before the start of the Windows operating system.
 - After Scandisk has been finished successfully Windows and the IOLMaster software will be started. Continue now as normal.

Rarely, switching off the device at the power switch may result in data loss or destruction of program modules.

In such a case, contact Carl Zeiss Meditec Service.

Checking the measurement functions

The test eyes (Fig. 73 und Fig. 74) supplied as accessories to the instrument serve for checking the proper function and calibration of the instrument.

Using these test eyes you can perform measurements as on a human eye. These measurements also require previous entry of name, first name and date of birth!

For checking the WTW measurement function (optional), use the provided calibration scale. Prior to the measurement, stick the scale onto the surface opposite the test eye (ACD) (applies only to later updates).


Advice

It is advisable to check the instrument calibration on a daily basis using the test eyes. You may print out the calibration results for documentation purposes.

If the measured values obtained with the test eyes are outside the specified tolerance range, please call Carl Zeiss Meditec or Authorized Service.

Remove holding pins for paper pads (7) beside the chin rest if inserted. Plug asymmetrical holder (6) into the corresponding holes. The test eye holder (1) is pivoted on pin (5) and may be fixed by means of fastening screw (4).

- On the **Options** menu, select *Test Eye*.
On the display, the word "Test Eye" appears.

- 1 Test eye holder
- 2 Test eye for AL and KER, labeled AL and R
- 3 Nominal values and tolerances
- 4 Fastening screw
- 5 Pin
- 6 Asymmetrical holder
- 7 Holding pins for paper pads (2x)
- 8 Chin rest

Fig. 73 Mounting the test eye

Axial length mode and keratometer

The test eye (2), labeled with AL, R and the corresponding nominal values and tolerances (3), serves for testing the axial length mode (AL) and the keratometer (R).


If the obtained measured values (with keratometer: the radius) are within the tolerance range (3) specified on the holder, the instrument is correctly calibrated.

Anterior chamber depth mode

The (larger) test eye (2) located at the side of the test eye holder (1), which is labeled with VKT (ACD), nominal value and tolerance, serves for testing the anterior chamber depth mode.

The structure of this lens simulates the cornea. Therefore, it must be clean and free of film. If necessary, wipe it using a dry cloth.

If the obtained measured values are within the specified tolerance range, the anterior chamber depth meter works correctly.

- 1 Test eye holder
- 2 Test eye for ACD; labeled VKT
- 3 Nominal values and tolerances
- 4 Asymmetrical holder
- 5 Mounting hole

Fig. 74 Testing the anterior-chamber depth mode

Explanation

When you test the axial length mode and the keratometer, it does not matter which side (right or left on simulated eye) is used, as the optical paths for the measurement are axisymmetric. When testing the anterior chamber depth meter, however, it is advisable to turn the asymmetrical holder (4) and perform the test also for the other side. For the comparison of right-hand and left-hand side, however, make sure that in both tests the test eye (2) is exactly perpendicular in front of the instrument.

View the video image to check if the adjustment criteria for an optimal optical section are met (--> *How to adjust the instrument, p. 71*).

Note:

The test eyes are excellently suitable also for training the operation of the IOLMaster.

The "Test Eye" mode can be deactivated again via **Options** menu, *Test eye*. If you enter a new patient (<N> or appropriate icon), this will be done automatically.

Checking the WTW function (optional)


The WTW scale (optional) (2, Fig. 75) serves for checking the precision of White-To-White measurements. The scale is to be stuck to the test eye holder in such a way, that its lines are vertical.

- Via the **Options** menu – *Test Eye*, select the "Test Eye" status and the WTW mode.
- Confirm the displayed instructions (Fig. 76) with **OK**
- Start a measurement:
If the measured values are within the tolerance range, the WTW function works correctly.

Note:

The WTW scale must fully fill the video window.
The scale (black lines) must be imaged sharply.

- 1 Test eye holder
- 2 WTW scale

Fig. 75 WTW scale (optional)


Fig. 76 Window instructing on the use of the WTW calibration scale

Remedy of printer troubles

Please, use printers recommended by Carl Zeiss Meditec only. The currently recommended printers are recallable on the Internet under:

http://www.meditec.zeiss.com/iol_master

If you encounter any problems during printing, first remove all printer drivers you do not use:

- On the **OPTIONS** menu, select *Setup, Add Printer*.


Fig. 77 Printer folder

- From the printers window select the icon for the connected printer (e. g. Epson Stylus Color 880). Right click to view the drop down menu. Select *Set as default*.
- Select the printers not connected (except for the "Add Printer" icon) and press the key. Follow the instructions displayed on the screen.
- Close the Printer folder again, when the uninstall routine is finished.

If the printer troubles should persist, please call our service representative.

In the United States call Customer Service Carl Zeiss Ophthalmic Division at 877-486-7473.

Instrument care

Warning:

Prior to cleaning turn off the instrument and disconnect the power cable!

Paper pads should be changed and chin rest and forehead rest should be disinfected after every patient.

For disinfection of parts with which the patient comes into contact (forehead and chin rest), use a disinfectant authorized for this application (see list of disinfectants and methods tested and approved by Robert-Koch-Institut, 14th edition, of 31 May, 2002).

You can clean all paint surfaces of the equipment with a moist, but not dripping cloth.

Do not use any aggressive or abrasive cleaning agents.

Wipe off sticky residues using a mild detergent and water.

To clean the display and the keyboard, use commercial cleaning cloths for computers and monitors.

Caution:

To avoid damage, absolutely take care that no moisture gets into the instrument or the keyboard.

Safety inspections

To ensure proper function of the unit, safety inspection of the unit should be made regularly. This technical safety inspection must be performed by authorized specialists once a year.

Consider the applicable local safety requirements.

Line voltage	100 ... 240 V AC ($\pm 10\%$)	
Frequency	50 ... 60 Hz	
Power consumption	90 VA	
Power connection	Power cord with power plug with earthing contact	
Enclosure Protection	IP 20	
Protection Class	I	
Device type	B (according to EN 60601-1)	
Dimensions		
Bench space	390 mm x 300 mm	
Height	max. 610 mm (headrest)	
Mass	approx. 18 kg	
Ambient conditions for intended use	Temperature:	+ 10 ... + 35 °C
	Relative humidity:	30 ... 75 %
	Air pressure:	800 ... 1060 hPa
Storage conditions	Temperature:	-10 ... + 55 °C
	Relative humidity:	10 ... 95 %
	Air pressure:	700 ... 1060 hPa
Shipment conditions (in original packing)	Temperature:	-40 ... + 70 °C
	Relative humidity:	10 ... 95 %
	Air pressure:	500 ... 1060 hPa
Measuring range		
Axial length		
Range	14 ... 40 mm	
Resolution of result display	0.01 mm	
Corneal radius		
Range	5 ... 10 mm (33 ... 67 D)	
Resolution of result display	0.01 mm	
Anterior chamber depth		
Range	1.5 ... 6.5 mm	
Resolution of result display	0.01 mm	
White To White (optional)		
Range	8 ... 16 mm	
Resolution of result display	0.1 mm	

Measurement information


Comparison results of IOLMaster versus conventional measurements at the human eye*:		
	Mean values of Differences	Standard Deviations
Axial length**	-0.01 mm	± 0.19 mm
Corneal radius***	-0.01 mm	± 0.05 mm
Anterior chamber depth**	0.07 mm	± 0.18 mm
IOLMaster measuring reproducibility*:		
	Regarding the human eye standard deviation****	regarding the artificial model eye standard deviation
Axial length	± 0.03 mm	± 0.01 mm
Corneal radius	± 0.02 mm	± 0.02 mm
Anterior chamber depth	± 0.10 mm	± 0.10 mm
* Data on File ** Compared to high precision ultrasound immersion device *** Compared to a handheld keratometer **** Standard deviation determined on basis of "x _i " calculated single standard deviations		

Optical radiation

Surrounding field illumination/ WTW determination	
Wavelength	880 nm
Delivered power	< 100 μ W
Axial length measurement	
Light source	Semiconductor diode laser (MMLD)
Wavelength	780 nm
Max. power for measurement	450 μ W
Max. power for alignment	80 μ W
Measuring time for individual measurement, Pulse width	0.5 s
Number of possible individual measurements	20 per eye and day
Laser class	
• In intended use	1 (EN 60825-1)
• Embedded Laser Class (not accessible)	3B
Fixation light for keratometer and anterior chamber depth measurement and WTW determination	
Wavelength	590 nm
Delivered power	< 1 μ W
Illumination for keratometer measurement	
Wavelength	880 nm
Delivered power	< 50 μ W
Slit illumination for anterior chamber depth measurement, integral irradiance	
UV (300 ... 400 nm)	0.00087 mW cm ⁻²
IR (700 ... 1100 nm)	0.04 mW cm ⁻²
	(in spectral region from 860 to 1,100 nm no emission of light source detectable)
L_B (phakic eye)	122.8 W (m ² sr) ⁻¹
L_A (aphakic eye)	125.5 W (m ² sr) ⁻¹

Spectral irradiance


spectral irradiance (white LED) in optical axis (0°) and in 33° -direction


"Optical axis" or "0°" corresponds to the direct view into the illuminating projector;
33° is the angle for the intended use in anterior-chamber depth measurement.

Spectrally assessed photochemical radiation densities

photochemical radiation density of slit illumination (white LED) for phacic eyes
 $L_B = 122,8 \text{ W/(m}^2 \text{ sr)}$

Photochemical radiation density of slit illumination (white LED) for aphakic eyes
 $L_A = 125,5 \text{ W/(m}^2 \text{ sr)}$

Note:

The spectrally assessed photochemical radiation densities L_B and L_A are a measure of the possibility of photochemical damage of the retina through light. L_B represents the measure for the phakic eye, L_A represents the measure for the aphakic eye or for the eyes of very young children. Values of L_B and L_A above $800 \text{ W/(m}^2 \text{ sr)}^{-1}$ are regarded as high. The radiation dose of the retina for a photochemical risk is calculated as the product of radiation density and exposure time. The recommended radiation dose is based on calculations of the American Conference of Governmental and Industrial Hygienists (ACGIH) Threshold Limit Values for Chemical Substances and Physical Agents (Edition: 1995-1996).

The measured photometric values of the IOLMaster are far below the levels that are regarded as high. Thus, the risk of damage through optical radiation is extremely low. Nevertheless, anterior chamber depth measurement with the IOLMaster should be limited to the time absolutely necessary for the diagnosis. The risk of damage may be higher, if fundus photography of the patient to be examined has been taken within the last 24 hours.

Technical details and delivery package subject to change.

Medical Product Act

The device meets the requirements of the EC Medical Device Directive 93/42/EEC and of their national equivalent in form of the German Medical Product Act (MPA).

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The screenshot shows the IOLMaster software interface. On the left, a list of patient records is displayed, including names like 'Muster, Max' and 'Mustermann, Doris'. The main area contains fields for patient information (Nom, Prénom, Date de Naissance, Date de l'examen, Numéro ID) and surgical details (Oeil opéré, AL, K1, K2, Opt. ACD, Impt. IOL (D), Réf. Post. Op., Date chirurgie, Date Post. Op.). Below these, there are sections for 'Base' and 'Nouveau' with various optical constants (Const. A, a0, pACD, SF) and their optimized values. A table on the right shows 'Données enregistrées' with values for different axial lengths (14..22mm, 22..25mm, 25..40mm) and the average AL. At the bottom, there are buttons for 'Charger', 'Nouveau', 'Optimiser', 'Effacer', 'OK', and 'Annuler'.

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The screenshot shows the 'Calcular IOL' (Calculate IOL) software window. It has a menu bar with 'Haigis', 'SRK II', 'HofnerQ', 'Holladay', 'SRK/T', and 'Cirugía refractiva anterior'. The patient name 'Mustermann, Felix, 02.05.1909' is at the top. The interface is split into two main columns for the right eye ('OD (derecha)') and the left eye ('OS (izquierda)'). Each column contains three sections: 'Valores medidos' (Measured values), 'Método de historia clínica' (Clinical history method), and 'Método de lente de contacto (CL)' (Contact lens method). The 'Valores medidos' section includes 'Longitud de eje' (axial length) and 'KER' (keratometry). The 'Método de historia clínica' section includes 'Poder refringente corneal preop.' (preoperative corneal refractive power), 'Vértice' (vertex distance), and pre/postoperative spherical and cylindrical refractive errors. The 'Método de lente de contacto' section includes 'Potencia CL' (CL power), 'Curva de base CL' (base curve CL), 'Refr. con CL' (refraction with CL), 'Refr. sin CL' (refraction without CL), and the calculated 'Poder refringente corneal' (corneal refractive power) and 'Radio' (radius). Checkboxes for 'Aplicar' (Apply) are present for the clinical history and contact lens methods. At the bottom are 'OK' and 'Abbrechen' (Cancel) buttons.

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Carl Zeiss Meditec AG